



# GOHOME™

## Patient Health Monitor



## GoHome™ — A telehealth communications hub for remote monitoring, communication, and care of patients at home

*The secure and simple way to connect patients to healthcare professionals.*

- Seamless bi-directional communication with the MTI Telehealth CarePortal™
- Automated reminders for medications, physiologic testing, and appointments
- Simplifies patient's diagnostic testing with an auto-pairing Bluetooth® connection to peripheral devices, and provides easy-to-follow testing instructions
- Automatically transfers physiologic data to healthcare providers via the MTI CarePortal
- Allows patients to connect and communicate with physicians, pharmacists, caregivers, and other health professionals, via secure video conferencing, email, and SMS messaging

**MONITOR. CONNECT. CARE.**



**MONITORED**  
THERAPEUTICS INC

## Optimize Patient Care and Outcomes



The *GoHome* connects patients with their physicians and caregivers, as well as with their families. Helping to maintain their health with remote monitoring and diagnostic testing enabled by automated reminders and easy-to-use communication tools.


The *GoHome* also gives patient families and caregivers the social tools to remain connected and involved from any distance.

## Automated Patient Reminders

Automated measurement, medication, and appointment reminders to improve treatment adherence. Medication Reminders include drug name, dose, time, and instructions.

Measurement and Medication Adherence Reports are sent to physicians or family members for remote tracking.

Patients can also receive and respond to subjective assessments through question sets and electronic patient diaries.



### For more information:

Email: [info@MonitoredRX.com](mailto:info@MonitoredRX.com)  
Visit: [www.MonitoredRX.com](http://www.MonitoredRX.com)

5940 Venture Dr., Suite C  
Dublin, OH 43017

3361-MTI-Rev A

## Simplified Patient Physiologic Data Collection

Simplifies patient's diagnostic testing with an auto-pairing Bluetooth® connection to peripheral devices.

The *GoHome* provides easy-to-follow, on-screen testing instructions for each type of device.

The *GoHome* collects physiologic data from:

- *GoSpiro*® Spirometer
- Blood Pressure Cuffs
- Weight Scales
- Pulse Oximeters
- Glucometers
- Other Physiologic Monitors

Once the test is completed, the *GoHome* automatically transfers the physiologic data using a HIPPA compliant protocol to healthcare providers via the MTI *CarePortal*.

## Easy-to-Use Technology

The *GoHome* comes with all software pre-loaded, all the patient needs to do is plug it in, and turn it on.

The *GoHome* is designed to look like a digital picture frame, not a medical device, ensuring that the unit will be kept ON (unlike computers and mobile devices).

## Easily Connect and Communicate

The *GoHome* enables easy and secure bi-directional communication with patients via email or SMS. Secure Video Chat with patients at the press of a single button.

The MTI *CarePortal* allows for easy communication with the *GoHome*, which is perfect for healthcare providers, caregivers or family members to monitor and assist the patient. Send simplified email communications with unique one-touch responses using CareText.™

## The Ideal Platform for:

- Hospitals and Health Systems
- Senior Care Service Companies
- Physician Group Practices
- Assisted Living Facilities
- Home Health Providers
- Clinical Trials
- Accountable Care Organizations (ACOs)

**MONITOR. CONNECT. CARE.**

